



OPTICS TECHNOLOGY

SINCE 1978

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2013 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

WILLEY MILL CE

INTERMEDIATE MODEL "ARTHUR THOMAS" TYPE



OT-WM-223

Electrically operated to work on 220 volts A.C Recommended for grinding of dry plant samples.

Also for a great variety of commercial dry materials. Unit having two stationary steel blades and rotary with steel edges provided in the cast iron chamber.

Willey mill is fitted on substantial cast iron base.

Model No.	WORKING CHAMBER
-----------	-----------------

OT-WM-223-A	40 mm x 25 mm
-------------	---------------

OT-WM-223-B	65 mm x 25 mm
-------------	---------------

OT-WM-223-C	75 mm x 25 mm
-------------	---------------

OT-WM-223-D	100 mm x 50 mm
-------------	----------------

OT-WM-223-E	200 mm x 75 mm
-------------	----------------

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in