



# OPTICS TECHNOLOGY

SINCE 1978

**Manufacturer & Exporters of : LAB TESTING INSTRUMENTS**



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010  
An ISO 10002:2018 CERTIFIED BS EN 12469:2000  
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

## AUTOMATIC TABLET HARDNESS TESTER

MICROPROCESSOR BASED



Digital tablet hardness tester are instruments for accuracy and detecting braking hardness of tablet. Simple and direct constructions made its outward appearance pleasing to the eyes.

### Working:

1. The tablet is loaded and pressed artificially.
2. It adopt high precision press sensor to ensure high test accuracy and repeatability.
3. It tests Tablet Hardness continuously. The biggest number is 100 pieces of a gang.
4. Time of testing will be saved by directly measuring the diameter.
5. Test data is display and locked automatically, so can test be reset and cycled.

### Specification:

<b>Hardness range</b>	: 2-199.9 N
<b>Hardness accuracy</b>	: Max. $\pm 0.5\%$
<b>Probe processing range</b>	: (3- 40)mm
<b>Maximum No.</b>	: In each group 100P
<b>Power Supply</b>	: 220V/ 50Hz/ 60W or 110V/ 60Hz/ 60W

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: [info@opticstechnology.in](mailto:info@opticstechnology.in)

[opticstechnology@gmail.com](mailto:opticstechnology@gmail.com)

Website : [www.opticstechnology.in](http://www.opticstechnology.in)