



SINCE 1978

OPTICS TECHNOLOGY

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022



BLOOD TUBE SHAKER (BLOOD ROLLER MIXER) CE



OT-BTS-158

Roller mixer is used for mixing blood samples, for prevention of blood coagulation and immune precipitation etc.

Features:

- * Gentle and highly efficient rocking and rolling motion.
- * Unique shape design
- * Four rollers with relevant space
- * Adjustable speed
- * Easy Maintenance

Technical Parameters:

Mixing Mode	: Rocking & Rolling
Speed Range	: 30-36RPM
Number of Roller	: 4pcs
Roller Length	: 230mm
Mixing Time	: 15-20 min
Power Supply	: 220V, 50/60Hz

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in