



OPTICS TECHNOLOGY

SINCE 1978

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

BREAK LOAD TESTER

CE



OT-BLT-341

Break Load Tester for lipstick is designed to known value of maximum load a lipstick can with stand before it break. The break load tester comes with display.

Technical Data:

Capacity	: 5 kg
Least count	: 0.001 kg
Motor	: 1/4 HP, single phase induction
Gear box	: Worm reduction type

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in