



SINCE 1978

OPTICS TECHNOLOGY

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

DEEP FREEZER (VERTICAL) CE



Temp. Range: (- 80°C)

- Puf insulated to prevent thermal losses
- Eco friendly refrigerant CFC free Heavy Duty Compressor
- Temp. controlled by Microprocessor based Temp. Controller (Digital LED Display)
- Very low noise
- Lock & key arrangement
- Door is sealed with Silicon Gasket
- Heavy duty Castor Wheels Mounted for easy movability
- Works on 220 VAC / 50Hz Mains.

Model No.	Capacity
OT-DF/V-135A	100 Ltrs.
OT-DF/V-135B	200 Ltrs.
OT-DF/V-135C	400 Ltrs.
OT-DF/V-135D	500 Ltrs.

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in