



OPTICS TECHNOLOGY

SINCE 1978

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

DESICCATOR



Glass Desiccator



Polycarbonate Desiccator

Plain Desiccators are physically the same as Vacuum Desiccator. Only difference is the absence of stop cock.
Grease should be applied on the flanges to make it air tight.

Material : Polycarbonate / Glass

DIA :

150 mm

210 mm

240 mm

300 mm

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in