



SINCE 1978

OPTICS TECHNOLOGY

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

GEL ROCKER SHAKER CE



Can be used for many applications like gel staining, culture aeration, bathing, blotting membranes etc. A metal platform with non-slip rubber mat moves in a gentle.

'See Saw' motion. The movement is controlled precisely with a DC motor. Cushion and rubber feet prevent unit from 'slipping'. Variable speed regulator

Motor : D.C.

Maximum Speed : 20 RPM

Speed controller : Analogue

Input Voltage : 230V 50Hz, A.C. Supply.

Platform Size :	OT-GRS-156-1	OT-GRS-156-2
	25 x 25 cms	30 x 30 cms

Optional

1. Digital RPM Indicator

2. Mechanical Timer

OR

Digital Timer

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in