



# OPTICS TECHNOLOGY

SINCE 1978

**Manufacturer & Exporters of : LAB TESTING INSTRUMENTS**



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010  
An ISO 10002:2018 CERTIFIED BS EN 12469:2000  
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022



## HORIZONTAL (RECIPROCATING) KAHN SHAKER CE



**OT-HKS-155**

D.C. Motor special system shaker to work on 220 Volts A.C. for variable and smooth function speed control by speed regulator. It is general purpose machine and found useful for mixing or agitation ingredients, solution in flasks, bottles, beakers and kahn racks etc. Upper tray covered with rubber floor and sides of the tray are slotted to take parallel and cross bars cover with rubber. The bars are tightened with washers and fly nuts which can be adjusted as per the requirement.

- It can carry any Flask of size 100ml, 250ml, 500ml.
- Body is made of mild steel duly powder coated

Model No.	Platform Size
<b>OT-HKS-155-A</b>	20" x 13" x 4"
<b>OT-HKS-155-B</b>	31" x 17" x 4"
<b>OT-HKS-155-C</b>	36" x 18" x 4"

### Optional:

1. Digital RPM Indicator
2. Mechanical Timer **OR** Digital Timer

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: [info@opticstechnology.in](mailto:info@opticstechnology.in)

[opticstechnology@gmail.com](mailto:opticstechnology@gmail.com)

Website : [www.opticstechnology.in](http://www.opticstechnology.in)