



SINCE 1978

OPTICS TECHNOLOGY



Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

HOT PLATE CE

(ROUND)



M.S. Body Cast Iron Top
with Energy Regulator
OT-HP-39 A (ER)



M.S. Body with S.S. Top
with Energy Regulator
OT-HP-39 B (ER)



S.S. Body with S.S. Top
with Energy Regulator
OT-HP-39 C (ER)

8" Dia Round

- Max. Temp. upto 350°C
- With Energy Regulator
- Plate has insulation
- Elements to work on 220/230 Volts A.C.

MOC:

- (a) M.S. Body with Cast Iron Top
- (b) M.S. Body with S.S. Top 316
- (c) S.S. Body with S.S. Top 316

Optional:

- Microprocessor Based PID Controller



M.S. Body Cast Iron Top
with PID Controller
OT-HP-39 D (PID)



M.S. Body with S.S. Top
with Energy Regulator
OT-HP-39 E (PID)



S.S. Body with S.S. Top
with PID Controller
OT-HP-39 F (PID)

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)
☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in,

opticstechnology@gmail.com,

Website : www.opticstechnology.in