



SINCE 1978

OPTICS TECHNOLOGY

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022



LABORATORY JACK



OT-LJ-274-MS

M.S. BODY



OT-LJ-274-SS

S.S. BODY

Lab Jacks are designed to support and adjust the height of a variety of laboratory apparatus such as hot plates, heating mantles, baths, and glassware, not for lifting loads.

Model No.	Size
OT-LJ-274 MS/SS-A	150 x 200 mm
OT-LJ-274 MS/SS-B	200 x 200 mm
OT-LJ-274 MS/SS-C	250 x 250 mm
OT-LJ-274 MS/SS-D	300 x 300 mm

MOC: MS-Mild steel duly powder coated SS-stainless steel
kindly select the one required in model number.

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in