



OPTICS TECHNOLOGY

SINCE 1978

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

MONOCULAR MICROSCOPE



OT-MMS-159

Monocular microscopes are often used in classrooms and laboratories for observing slide samples and it is an instrument for viewing small things through a single lens.

Stereo microscopes are available as monocular microscopes

- Illumination by Plano concave mirror in fork mount.
- Objectives : 10x, 45x
- Optical combination
- Eye Piece : 10x, 40x (Huygenian)

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in