



OPTICS TECHNOLOGY

SINCE 1978

Manufacturer & Exporters of : LAB TESTING INSTRUMENTS



An ISO 9001:2015 CERTIFIED IEC 61010-1:2010
An ISO 10002:2018 CERTIFIED BS EN 12469:2000
An ISO 13485:2016 CERTIFIED ISO 29463-5:2022

PARAFIN WAX BATH

CE



OT-PWB-240

Double walled outer made of M.S. Duly Powder Coated / S.S. 304

Inner stainless steel 304 grade.

The bath consist of 2 pilot lamps & on-off switch.

Supplied with immersion heaters.

Temp. controlled by Thermostat from ambient temp.+5°C to 95°C $\pm$ 1°C

Model No.	Capacity
OT-PWB-240-MS/SS-A	2 Ltrs.
OT-PWB-240-MS/SS-B	3 Ltrs.
OT-PWB-240-MS/SS-C	4 Ltrs.
OT-PWB-240-MS/SS-D	5 Ltrs.

2/3, DSIDC, (Behind Fun Cinema) NSP, Pitampura, Delhi-110034 (INDIA)

☎ 9811238150, 9811238151, 9811236960

E-mail: info@opticstechnology.in

opticstechnology@gmail.com

Website : www.opticstechnology.in