



ROTARY FLASKS SHAKER (MECHANICAL SHAKER) CE



Rotates specimens in horizontal plane in 3/4" circles, constructed of mild steel finished with stoved enamel point, Turned up edges and rubber sheeting on platform to prevent slipping of specimen containers, platform rigidly supported by steel legs for smooth and silent rotation, mounted on four rubber suction feet to prevent creeping.

Models No.

Size

OT-RFS-153/A	12" x 12"	(9 Flask of 250 ml)	OR Mix
OT-RFS-153/B	18" x 18"	(16 Flask of 250 ml)	OR Mix
OT-RFS-153/C	24" x 24"	(25 Flask of 250 ml)	OR Mix
OT-RFS-153/D	36" x 36"	(36 Flask of 250 ml)	OR Mix
OT-RFS-153/E	48" x 48"	(48 Flask of 250 ml)	OR Mix
OT-RFS-153/F	60" x 60"	(60 Flask of 250 ml)	OR Mix

OR Any Non Standard size as per requirement

Optional :

1. Digital R.P.M. indicator
2. Digital Timer

MOC :

- (a) M.S. Top with Rubber Disc
- (b) M.S. Top with S.S. Lotus Clamps
- (c) S.S. Top with S.S. Lotus Clamps